

Case Competition			
Registration No.	Name	Affiliation	Case Title
C010	Junpei Koike	Kyoto Katsura Hospital, Japan	A case of chronic total occlusion of the right coronary artery in which the wire stuck to the microcatheter after rendezvous
C011	Koki Omi	Nihonkai General Hospital, Japan	Bailout case for the pitfall of Knuckle Wire Technique (KWT) in long CTO PCI
C014	Alexander Ijsselmuiden	Amphia Hospital, The Netherlands	Bail out treatment through a CTO
C018	Ali Raza	University of Alberta, Canada	A Novel Technique for Tackling Trapped Devices
Poster Presentation			
Registration No.	Name	Affiliation	Case Title
C001	Jae Young Cho	Wonkwang University Hospital, Korea	CTO Intervention of Distal RCA bifurcation with Poor Collateral Flow
C002	Kiwamu Sudo	Saiseikai Utsunomiya Hospital, Japan	A case of successful bailout for coronary perforation during antegrade wiring
C003	Takuma Tsuda	Nagoya Ekisaikai Hospital, Japan	Stentless PCI for LAD CTO by DCA
C004	Yuya Nakagawa	Hakodate Medical Association Hospital, Japan	Complications when using snare catheter to catch retrograde wire
C005	Keisuke Nakabayashi	Kasukabe Chuo General Hospital, Japan	The wire rendezvous and chasing wire technique in the bidirectional approach for the percutaneous coronary intervention for chronic total occlusion with a single guiding catheter
C006	Marohito Nakata	Naha City Hospital, Japan	Success bails out of A-C channel perforation by coil embolization
C007	Kook-Jin Chun	Pusan National University Yangsan Hospital, Korea	Antegrade Approach to LAD os CTO Jailed by Previously Deployed Stent
C008	Jin Sup Park	Pusan National University Hospital, Korea	Nightmare of the previous implanted stent for treating CTO segment
C009	Kohei Asada	Shiga University of Medical Science, Japan	A Nightmare of LAD CTO with Ellis Type III Coronary Perforation
C012	Jong-Hwa Ahn	Gyeongsang National University Changwon Hospital, Korea	Failed reverse-CART technique case: that is because retrograde microcatheter was not crossing
C013	Yuta Imai	Saiseikai Shiga Prefecture Hospital, Japan	Entrapment of Gaia Next in Heavy Calcified RCA CTO
C015	Takuya Sakemi	Iizuka Hospital, Japan	A pitfall of rCART in a LAD CTO case
C016	Min Gyu Kang	Gyeongsang National University Hospital, Korea	Coronary perforation complicating cardiac tamponade during complex percutaneous coronary intervention
C017	Munenori Okubo	Gifu Heart Center, Japan	My worst case; Cardiac tamponade 2days after CTO PCI
C020	Motosu Ando	Okamura Memorial Hospital, Japan	A pitfall just before the end
C021	Makoto Muto	Saitama Prefectural Cardiovascular and Respiratory Center, Japan	Guide wire entrapment in PCI to RCA ostial ISR CTO
C022	Hirohisa Endo	Juntendo University, Graduate School of Medicine, Japan	Successful retrieval of the dislodged ultrathin strut stent (Orsiro stent) using a gooseneck snare kit
C023	Yuki Suzuka	Mimihara General Hospital, Japan	A case of RCA-CTO with high bend and calcification, complicated with perforation of epicardial channel
C024	Ritsu Yoshida	Kyoto Katsura Hospital, Japan	A case of difficult to determine the cause of stuck Tornus in RCA CTO
C025	Kenichi Oota	Tachibanadai Hospital, Japan	A LAD-CTO case that dissection occurred in the LMT

The 20th CTO Club Poster Posting

1. Posting & Removal

Location: Poster Posting Area in Exhibition Hall (5F, WINC AICHI)


Posting: Friday, June 14, 8:00-12:00

Removal: Saturday, June 15, 13:30-16:00

- Please post and remove your posters during the periods shown above. If you cannot come during the posting period, please post your poster as soon as you arrive at the venue.
- Posters which are left on the panel after the removal time will be removed and kept at the CTO Club administration office, and it will be disposed after the meeting.

2. Poster Posting

- Size of panel: W90cm, H210cm
- Poster numbers are prepared by the CTO Club
- Presenters are requested to prepare the abstract title, institution, presenter's name, co-authors' institution and co-authors' name
- Push pins will be provided at the Poster Posting Area
- Posters must be written in English



第 20 回 CTO Club ポスター掲示要項

1. ポスター貼付及び撤去時間

場所： 展示会場内 Poster Posting Area（ウインクあいち 5F）

貼付： 6 月 14 日（金） 8:00～12:00

撤去： 6 月 15 日（土） 13:30～16:00

- 上記時間内に各自でポスター貼付及び撤去をお願いいたします。
貼付時間内での貼付が難しい場合は、会場に到着次第、貼付をお願いいたします。
- 撤去時間を過ぎても掲示されているポスターにつきましては、事務局にて撤去、保管いたします。会期終了後には事務局にて処分いたしますので予めご了承下さい。

2. ポスター掲示について

- パネルの大きさは幅 90cm、高さ 210cm です。
- 演題番号は事務局で準備いたします。
- 演題、所属、発表者名は各自で準備してください。
- ポスターの貼付に必要な画鋏はパネルに用意しています。
- ポスターの言語は英語です。

